




EDNRA Antibody

Product Code	CSB-PA002447
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25101
Immunogen	Synthesized peptide derived from the C-terminal region of Human ETAR.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Monkey
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	EDNRA; ETA; ETRA; Endothelin-1 receptor; Endothelin A receptor; ET-A; ETA-R; hET-AR
Immunogen Species	Homo sapiens (Human)
Target Names	EDNRA

Image



Western Blot analysis of COS7 cells using ETAR Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.