





EZR Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA002456 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P15311 |
| Immunogen | Synthesized peptide derived from Human Ezrin around the non-phosphorylation site of T566. |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | WB, ELISA; WB:1:500-1:2000, ELISA:1:20000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | EZR; VIL2; Ezrin; Cytovillin; Villin-2; p81 |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | EZR |

Image



Western Blot analysis of HepG2 cells using Ezrin Polyclonal Antibody



Western Blot analysis of Jurkat cells using Ezrin Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.