




FADD Antibody

Product Code	CSB-PA002463
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13158
Immunogen	Synthesized peptide derived from Human FADD around the non-phosphorylation site of S194.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	FADD; MORT1; GIG3; Protein FADD; FAS-associated death domain protein; FAS-associating death domain-containing protein; Growth-inhibiting gene 3 protein; Mediator of receptor induced toxicity
Immunogen Species	Homo sapiens (Human)
Target Names	FADD

Image



Western Blot analysis of HeLa cells using FADD Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.