




RIMKLA Antibody

Product Code	CSB-PA002477
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IXN7
Immunogen	Synthesized peptide derived from the Internal region of Human FAM80A.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Monkey
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	RIMKLA; FAM80A; N-acetylaspartyl-glutamate synthetase A; NAAG synthetase A; NAAGS; Ribosomal protein S6 modification-like protein A
Immunogen Species	Homo sapiens (Human)
Target Names	RIMKLA

Image



Western Blot analysis of Jurkat cells using FAM80A Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.