




FASTKD1 Antibody

Product Code	CSB-PA002493
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q53R41
Immunogen	Synthesized peptide derived from the Internal region of Human FASTKD1.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Monkey
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	FASTKD1; KIAA1800; FAST kinase domain-containing protein 1
Immunogen Species	Homo sapiens (Human)
Target Names	FASTKD1

Image



Western Blot analysis of COS7 cells using FASTKD1 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.