




FASN Antibody

Product Code	CSB-PA002507
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49327
Immunogen	Synthesized peptide derived from the Internal region of Human Fatty Acid Synthase.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	[Acyl-carrier-protein] S acetyltransferase antibody; [Acyl-carrier-protein] S malonyltransferase antibody; 3-hydroxypalmitoyl-[acyl-carrier-protein] dehydratase antibody; 3-oxoacyl-[acyl-carrier-protein] reductase antibody; 3-oxoacyl-[acyl-carrier-protein] synthase antibody; Enoyl-[acyl-carrier-protein] reductase antibody; FAS antibody; FAS_HUMAN antibody; FASN antibody; Fatty acid synthase antibody; MGC14367 antibody; MGC15706 antibody; OA 519 antibody; Oleoyl-[acyl-carrier-protein] hydrolase antibody; SDR27X1 antibody; Short chain dehydrogenase/reductase family 27X member 1 antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FASN

Image



Western Blot analysis of A549 cells using Fatty Acid Synthase Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.