




FASN Antibody

Product Code	CSB-PA002507
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49327
Immunogen	Synthesized peptide derived from the Internal region of Human Fatty Acid Synthase.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	[Acyl-carrier-protein] S acetyltransferase antibody; [Acyl-carrier-protein] S malonyltransferase antibody; 3-hydroxypalmitoyl-[acyl-carrier-protein] dehydratase antibody; 3-oxoacyl-[acyl-carrier-protein] reductase antibody; 3-oxoacyl-[acyl-carrier-protein] synthase antibody; Enoyl-[acyl-carrier-protein] reductase antibody; FAS antibody; FAS_HUMAN antibody; FASN antibody; Fatty acid synthase antibody; MGC14367 antibody; MGC15706 antibody; OA 519 antibody; Oleoyl-[acyl-carrier-protein] hydrolase antibody; SDR27X1 antibody; Short chain dehydrogenase/reductase family 27X member 1 antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FASN

Image



Western Blot analysis of A549 cells using Fatty Acid Synthase Polyclonal Antibody