





FER Antibody

Product Code	CSB-PA002517
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16591
Immunogen	Synthesized peptide derived from Human Fer around the non-phosphorylation site of Y402.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	FER; TYK3; Tyrosine-protein kinase Fer; Feline encephalitis virus-related kinase FER; Fujinami poultry sarcoma/Feline sarcoma-related protein Fer; Proto-oncogene c-Fer; Tyrosine kinase 3; p94-Fer
Immunogen Species	Homo sapiens (Human)
Target Names	FER

Image



Western Blot analysis of HeLa cells using Fer Polyclonal Antibody



Western Blot analysis of NIH-3T3 cells using Fer Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.