




FGL2 Antibody

Product Code	CSB-PA002521
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14314
Immunogen	Synthesized peptide derived from the N-terminal region of Human Fgl2.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	FGL2; Fibroleukin; Fibrinogen-like protein 2; pT49
Immunogen Species	Homo sapiens (Human)
Target Names	FGL2

Image



Western Blot analysis of 293 cells using Fgl2 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.