




FLNA Antibody

Product Code	CSB-PA002535
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21333
Immunogen	Synthesized peptide derived from Human Filamin 1 around the non-phosphorylation site of S2152.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	FLNA; FLN; FLN1; Filamin-A; FLN-A; Actin-binding protein 280; ABP-280; Alpha-filamin; Endothelial actin-binding protein; Filamin-1; Non-muscle filamin
Immunogen Species	Homo sapiens (Human)
Target Names	FLNA

Image



Western Blot analysis of HuvEc cells using Filamin 1 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.