





FN1 Antibody

Product Code	CSB-PA002552
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02751
Immunogen	Synthesized peptide derived from the C-terminal region of Human FN1.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	FN1; FN; Fibronectin; FN; Cold-insoluble globulin; CIG
Immunogen Species	Homo sapiens (Human)
Target Names	FN1


Image



Western Blot analysis of HeLa cells using FN1 Polyclonal Antibody



Western Blot analysis of 293T cells using FN1 Polyclonal Antibody



Western blot analysis of 3T3 KB K562 293T lysis using FN1 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.