




PLK3 Antibody

Product Code	CSB-PA002553
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H4B4
Immunogen	Synthesized peptide derived from the Internal region of Human Fnk.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PLK3; CNK; FNK; PRK; Serine/threonine-protein kinase PLK3; Cytokine-inducible serine/threonine-protein kinase; FGF-inducible kinase; Polo-like kinase 3; PLK-3; Proliferation-related kinase
Immunogen Species	Homo sapiens (Human)
Target Names	PLK3

Image



Western Blot analysis of HepG2 cells using Fnk Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.