




FOXD3 Antibody

Product Code	CSB-PA002566
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UJU5
Immunogen	Synthesized peptide derived from the Internal region of Human FoxD3.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	FOXD3; HFH2; Forkhead box protein D3; HNF3/FH transcription factor genesis
Immunogen Species	Homo sapiens (Human)
Target Names	FOXD3

Image



Western Blot analysis of 293 cells using FoxD3 Polyclonal Antibody