





FOXO1 Antibody

Product Code	CSB-PA002569
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12778
Immunogen	Synthesized peptide derived from Human FoxO1 around the non-phosphorylation site of S319.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	FOXO1; FKHR; FOXO1A; Forkhead box protein O1; Forkhead box protein O1A; Forkhead in rhabdomyosarcoma
Immunogen Species	Homo sapiens (Human)
Target Names	FOXO1

Image



Western Blot analysis of HeLa cells using FoxO1 Polyclonal Antibody



Western Blot analysis of NIH-3T3 RAT-MUSCLE cells using FoxO1 Polyclonal Antibody