





FOXO1 Antibody

Product Code	CSB-PA002571
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12778
Immunogen	Synthesized peptide derived from Human FoxO1A around the non-phosphorylation site of S329.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	FOXO1; FKHR; FOXO1A; Forkhead box protein O1; Forkhead box protein O1A; Forkhead in rhabdomyosarcoma
Immunogen Species	Homo sapiens (Human)
Target Names	FOXO1

Image



Western Blot analysis of HeLa cells using FoxO1A Polyclonal Antibody



Western Blot analysis of HELA Jurkat cells using FoxO1A Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.