





FYB Antibody

Product Code	CSB-PA002589
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15117
Immunogen	Synthesized peptide derived from the C-terminal region of Human Fyb.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	FYB; SLAP130; FYN-binding protein; Adhesion and degranulation promoting adaptor protein; ADAP; FYB-120/130; p120/p130; FYN-T-binding protein; SLAP-130; SLP-76-associated phosphoprotein
Immunogen Species	Homo sapiens (Human)
Target Names	FYB1

Image



Western Blot analysis of COLO205 cells using Fyb Polyclonal Antibody



Western Blot analysis of Jurkat cells using Fyb Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.