




G3BP1 Antibody

Product Code	CSB-PA002599
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13283
Immunogen	Synthesized peptide derived from Human G3BP1 around the non-phosphorylation site of S232.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Monkey
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	G3BP1; G3BP; Ras GTPase-activating protein-binding protein 1; G3BP-1; ATP-dependent DNA helicase VIII; hDH VIII; GAP SH3 domain-binding protein 1
Immunogen Species	Homo sapiens (Human)
Target Names	G3BP1

Image



Western Blot analysis of HeLa cells using G3BP1 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.