





GAB2 Antibody

Product Code	CSB-PA002610
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UQC2
Immunogen	Synthesized peptide derived from the Internal region of Human Gab 2.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GAB2; KIAA0571; GRB2-associated-binding protein 2; GRB2-associated binder 2; Growth factor receptor bound protein 2-associated protein 2; pp100
Immunogen Species	Homo sapiens (Human)
Target Names	GAB2

Image



Western Blot analysis of K562 cells using Gab 2 Polyclonal Antibody



Western Blot analysis of K562 cells using Gab 2 Polyclonal Antibody