





LGALS9 Antibody

Product Code	CSB-PA002634
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00182
Immunogen	Synthesized peptide derived from the Internal region of Human Galectin-9.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	LGALS9; Galectin-9; Gal-9; Ecalectin; Tumor antigen HOM-HD-21
Immunogen Species	Homo sapiens (Human)
Target Names	LGALS9

Image



Western Blot analysis of HepG2 cells using Galectin-9 Polyclonal Antibody



Western Blot analysis of VEC cells using Galectin-9 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.