




GAP43 Antibody

Product Code	CSB-PA002640
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17677
Immunogen	Synthesized peptide derived from Human GAP-43 around the non-phosphorylation site of S41.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, ELISA; IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GAP43; Neuromodulin; Axonal membrane protein GAP-43; Growth-associated protein 43; Neural phosphoprotein B-50; pp46
Immunogen Species	Homo sapiens (Human)
Target Names	GAP43

Image



Western Blot analysis of HELA KB cells using GAP-43 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.