




GATA1 Antibody

Product Code	CSB-PA002643
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15976
Immunogen	Synthesized peptide derived from Human GATA-1 around the non-phosphorylation site of S142.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IP, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IP:2-5ug, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GATA1; ERYF1; GF1; Erythroid transcription factor; Eryf1; GATA-binding factor 1; GATA-1; GF-1; NF-E1 DNA-binding protein
Immunogen Species	Homo sapiens (Human)
Target Names	GATA1

Image



Western Blot analysis of HT29 cells using GATA-1 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.