

🕜 Tel: +1-301-363-4651 🛛 🖂 Email: cusabio@cusabio.com 🥃 Website: www.cusabio.com 🍙

GATA1 Antibody

Product Code	CSB-PA002643	
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.	
Uniprot No.	P15976	
Immunogen	Synthesized peptide derived from Human GATA-1 around the non- phosphorylation site of S142.	
Raised In	Rabbit	
Species Reactivity	Human,Mouse,Rat	
Tested Applications	WB, IHC, IP, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IP:2-5ug, IF:1:200-1:1000, ELISA:1:10000	
Form	Liquid	
Conjugate	Non-conjugated	
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.	
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity- chromatography using epitope-specific immunogen.	
Isotype	IgG	
Alias	GATA1; ERYF1; GF1; Erythroid transcription factor; Eryf1; GATA-binding factor 1; GATA-1; GF-1; NF-E1 DNA-binding protein	
Immunogen Species	Homo sapiens (Human)	
Target Names	GATA1	
Image		m Blot analysis of HT29 cells using 1 Polyclonal Antibody

1