




GIMAP5 Antibody

Product Code	CSB-PA002686
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96F15
Immunogen	Synthesized peptide derived from the C-terminal region of Human GIMAP5.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GIMAP5; IAN4L1; IAN5; IMAP3; GTPase IMAP family member 5; Immunity-associated nucleotide 4-like 1 protein; Immunity-associated nucleotide 5 protein; IAN-5; hIAN5; Immunity-associated protein 3
Immunogen Species	Homo sapiens (Human)
Target Names	GIMAP5

Image



Western Blot analysis of HT29 cells using GIMAP5 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.