




GK Antibody

Product Code	CSB-PA002694
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P32189
Immunogen	Synthesized peptide derived from the C-terminal region of Human GK1.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GK; Glycerol kinase; GK; Glycerokinase; ATP:glycerol 3-phosphotransferase
Immunogen Species	Homo sapiens (Human)
Target Names	GK

Image



Western Blot analysis of Jurkat cells using GK1 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.