





GLP1R Antibody

Product Code	CSB-PA002699
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43220
Immunogen	Synthesized peptide derived from the Internal region of Human GLP-1R.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GLP1R; Glucagon-like peptide 1 receptor; GLP-1 receptor; GLP-1-R; GLP-1R
Immunogen Species	Homo sapiens (Human)
Target Names	GLP1R

Image



Western Blot analysis of HepG2 cells using GLP-1R Polyclonal Antibody



Western Blot analysis of HepG2 HuvEc cells using GLP-1R Polyclonal Antibody

