





GCGR Antibody

Product Code	CSB-PA002705
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P47871
Immunogen	Synthesized peptide derived from the Internal region of Human Glucagon Receptor.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GCGR; Glucagon receptor; GL-R
Immunogen Species	Homo sapiens (Human)
Target Names	GCGR

Image



Western Blot analysis of COLO205 cells using Glucagon Receptor Polyclonal Antibody



Western Blot analysis of COLO205 cells using
Glucagon Receptor Polyclonal Antibody