




B4GALNT1 Antibody

Product Code	CSB-PA002730
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q00973
Immunogen	Synthesized peptide derived from the Internal region of Human GM2/GD2 synthase.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	B4GALNT1; GALGT; SIAT2; Beta-1; 4 N-acetylgalactosaminyltransferase 1; (N-acetylneuraminy)-galactosylglucosylceramide; GM2/GD2 synthase; GalNAc-T
Immunogen Species	Homo sapiens (Human)
Target Names	B4GALNT1

Image



Western Blot analysis of HeLa cells using GM2/GD2 synthase Polyclonal Antibody