





GNL3L Antibody

Product Code	CSB-PA002732
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NVN8
Immunogen	Synthesized peptide derived from the Internal region of Human GNL3L.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GNL3L; Guanine nucleotide-binding protein-like 3-like protein
Immunogen Species	Homo sapiens (Human)
Target Names	GNL3L


Image



Western Blot analysis of HT29 cells using GNL3L Polyclonal Antibody



Western Blot analysis of NIH-3T3 LOVO cells using GNL3L Polyclonal Antibody



Western blot analysis of mouse-brain HELA KB SH-SY5Y 293T 3T3 Mouse-kidney lysis using GNL3L antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.