




GPR101 Antibody

Product Code	CSB-PA002738
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96P66
Immunogen	Synthesized peptide derived from the C-terminal region of Human GPR101.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Monkey
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GPR101; Probable G-protein coupled receptor 101
Immunogen Species	Homo sapiens (Human)
Target Names	GPR101

Image



Western Blot analysis of COS7 cells using GPR101 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.