




P2RY14 Antibody

Product Code	CSB-PA002743
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15391
Immunogen	Synthesized peptide derived from the Internal region of Human GPR105.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	P2RY14; GPR105; KIAA0001; P2Y purinoceptor 14; P2Y14; G-protein coupled receptor 105; UDP-glucose receptor
Immunogen Species	Homo sapiens (Human)
Target Names	P2RY14

Image



Western Blot analysis of HepG2 cells using GPR105 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.