





GPR110 Antibody

Product Code	CSB-PA002753
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5T601
Immunogen	Synthesized peptide derived from the C-terminal region of Human GPR110.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GPR110; PGR19; Probable G-protein coupled receptor 110; G-protein coupled receptor KPG_012; G-protein coupled receptor PGR19
Immunogen Species	Homo sapiens (Human)
Target Names	ADGRF1

Image



Western Blot analysis of 3T3 cells using GPR110 Polyclonal Antibody



Western Blot analysis of NIH-3T3 K562 cells using GPR110 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.