





# FFAR4 Antibody

<b>Product Code</b>	CSB-PA002756
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q5NUL3
<b>Immunogen</b>	Synthesized peptide derived from the Internal region of Human GPR120.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Isotype</b>	IgG
<b>Alias</b>	O3FAR1; GPR120; GPR129; PGR4; Omega-3 fatty acid receptor 1; G-protein coupled receptor 120; G-protein coupled receptor 129; G-protein coupled receptor GT01; G-protein coupled receptor PGR4
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	FFAR4

## Image



Western Blot analysis of HEPG2-UV cells using GPR120 Polyclonal Antibody



Western Blot analysis of HepG2-UV cells using GPR120 Polyclonal Antibody