





FFAR4 Antibody

Product Code	CSB-PA002756
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5NUL3
Immunogen	Synthesized peptide derived from the Internal region of Human GPR120.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	O3FAR1; GPR120; GPR129; PGR4; Omega-3 fatty acid receptor 1; G-protein coupled receptor 120; G-protein coupled receptor 129; G-protein coupled receptor GT01; G-protein coupled receptor PGR4
Immunogen Species	Homo sapiens (Human)
Target Names	FFAR4

Image



Western Blot analysis of HEPG2-UV cells using GPR120 Polyclonal Antibody



Western Blot analysis of HepG2-UV cells using GPR120 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.