





GPR152 Antibody

Product Code	CSB-PA002760
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TDT2
Immunogen	Synthesized peptide derived from the N-terminal region of Human GPR152.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GPR152; PGR5; Probable G-protein coupled receptor 152; G-protein coupled receptor PGR5
Immunogen Species	Homo sapiens (Human)
Target Names	GPR152

Image



Western Blot analysis of k562 cells using GPR152 Polyclonal Antibody



Western Blot analysis of various cells using GPR152 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.