




GPR171 Antibody

Product Code	CSB-PA002766
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14626
Immunogen	Synthesized peptide derived from the Internal region of Human GPR171.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GPR171; H963; Probable G-protein coupled receptor 171; G-protein coupled receptor H963
Immunogen Species	Homo sapiens (Human)
Target Names	GPR171

Image



Western Blot analysis of HT29 cells using GPR171 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.