




GPR83 Antibody

Product Code	CSB-PA002792
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NYM4
Immunogen	Synthesized peptide derived from the C-terminal region of Human GPR83.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GPR83; GPR72; KIAA1540; Probable G-protein coupled receptor 83; G-protein coupled receptor 72
Immunogen Species	Homo sapiens (Human)
Target Names	GPR83

Image



Western Blot analysis of HuvEc cells using GPR83 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.