




LPAR5 Antibody

Product Code	CSB-PA002795
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H1C0
Immunogen	Synthesized peptide derived from the Internal region of Human GPR92.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	LPAR5; GPR92; GPR93; Lysophosphatidic acid receptor 5; LPA receptor 5; LPA-5; G-protein coupled receptor 92; G-protein coupled receptor 93
Immunogen Species	Homo sapiens (Human)
Target Names	LPAR5

Image



Western Blot analysis of K562 cells using GPR92 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.