




GORASP2 Antibody

Product Code	CSB-PA002804
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H8Y8
Immunogen	Synthesized peptide derived from the Internal region of Human GRASP55.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GORASP2; GOLPH6; Golgi reassembly-stacking protein 2; GRS2; Golgi phosphoprotein 6; GOLPH6; Golgi reassembly-stacking protein of 55 kDa; GRASP55; p59
Immunogen Species	Homo sapiens (Human)
Target Names	GORASP2

Image



Western Blot analysis of LOVO cells using GRASP55 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.