





GPRIN1 Antibody

Product Code	CSB-PA002812
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z2K8
Immunogen	Synthesized peptide derived from the Internal region of Human GRIN1.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GPRIN1; KIAA1893; G protein-regulated inducer of neurite outgrowth 1; GRIN1
Immunogen Species	Homo sapiens (Human)
Target Names	GPRIN1

Image



Western Blot analysis of HT29 cells using GRIN1 Polyclonal Antibody



Western Blot analysis of PC3 KB cells using GRIN1 Polyclonal Antibody

