





ADRBK1 Antibody

Product Code	CSB-PA002823
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25098
Immunogen	Synthesized peptide derived from Human GRK 2 around the non-phosphorylation site of Y86.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ADRBK1; BARK; BARK1; GRK2; Beta-adrenergic receptor kinase 1; Beta-ARK-1; G-protein coupled receptor kinase 2
Immunogen Species	Homo sapiens (Human)
Target Names	GRK2

Image



Western Blot analysis of HuvEc cells using GRK 2 Polyclonal Antibody



Western Blot analysis of K562 HuvEc cells using GRK 2 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.