





# GRK5 Antibody

<b>Product Code</b>	CSB-PA002827
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P34947
<b>Immunogen</b>	Synthesized peptide derived from the Internal region of Human GRK 5.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Isotype</b>	IgG
<b>Alias</b>	GRK5; GPRK5; G protein-coupled receptor kinase 5; G protein-coupled receptor kinase GRK5
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	GRK5

## Image



Western Blot analysis of COLO205 cells using GRK 5 Polyclonal Antibody



Western Blot analysis of Jurkat K562 NIH-3T3 HELA HepG2 HuvEc COLO cells using GRK 5 Polyclonal Antibody

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.