





GSK3A Antibody

Product Code	CSB-PA002836
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49840
Immunogen	Synthesized peptide derived from Human GSK3 α around the non-phosphorylation site of S21.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, IP, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IP:2-5ug, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GSK3A; Glycogen synthase kinase-3 alpha; GSK-3 alpha; Serine/threonine-protein kinase GSK3A
Immunogen Species	Homo sapiens (Human)
Target Names	GSK3A

Image



Western Blot analysis of Jurkat cells using GSK3 α Polyclonal Antibody



Western Blot analysis of Jurkat cells using GSK3α Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.