





GSK3B Antibody

Product Code	CSB-PA002848
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49841
Immunogen	Synthesized peptide derived from Human GSK3 β around the non-phosphorylation site of S9.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, IP, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IP:2-5ug, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GSK3B; Glycogen synthase kinase-3 beta; GSK-3 beta; Serine/threonine-protein kinase GSK3B
Immunogen Species	Homo sapiens (Human)
Target Names	GSK3B


Image



Western Blot analysis of 293 cells using GSK3 β Polyclonal Antibody



Western Blot analysis of MOUSE-BRAIN cells using GSK3 β Polyclonal Antibody



Western Blot analysis of NIH-3T3 cells using GSK3 β Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.