




HCCS Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA002862 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P53701 |
| Immunogen | Synthesized peptide derived from the Internal region of Human HCCS. |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Monkey |
| Tested Applications | WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | HCCS; CCHL; Cytochrome c-type heme lyase; CCHL; Holocytochrome c-type synthase |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | HCCS |

Image



Western Blot analysis of COS7 cells using HCCS Polyclonal Antibody