





HDAC3 Antibody

Product Code	CSB-PA002880
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15379
Immunogen	Synthesized peptide derived from Human HDAC3 around the non-phosphorylation site of S424.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	HDAC3; Histone deacetylase 3; HD3; RPD3-2; SMAP45
Immunogen Species	Homo sapiens (Human)
Target Names	HDAC3

Image



Western Blot analysis of hela cells using HDAC3 Polyclonal Antibody



Western Blot analysis of HELA cells using HDAC3 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.