





HDAC5/HDAC9 Antibody

Product Code	CSB-PA002890
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UQL6/Q9UKV0
Immunogen	Synthesized peptide derived from Human HDAC5/9 around the non-phosphorylation site of S259/220.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	HDAC5; KIAA0600; Histone deacetylase 5; HD5; Antigen NY-CO-9; HDAC9; HDAC7; HDAC7B; HDRP; KIAA0744; MITR; Histone deacetylase 9; HD9; Histone deacetylase 7B; HD7; HD7b; Histone deacetylase-related protein; MEF2-interacting transcription rep
Immunogen Species	Homo sapiens (Human)
Target Names	HDAC5/HDAC9

Image



Western Blot analysis of HepG2 cells using HDAC5/9 Polyclonal Antibody



Western Blot analysis of NIH-3T3 RAT-intestines cells using HDAC5/9 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.