




HDAC6 Antibody

Product Code	CSB-PA002895
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBN7
Immunogen	Synthesized peptide derived from Human HDAC6 around the non-phosphorylation site of S22.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	HDAC6; KIAA0901; JM21; Histone deacetylase 6; HD6
Immunogen Species	Homo sapiens (Human)
Target Names	HDAC6

Image



Western Blot analysis of HepG2 cells using HDAC6 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.