





# HDAC7 Antibody

<b>Product Code</b>	CSB-PA002897
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q8WUI4
<b>Immunogen</b>	Synthesized peptide derived from Human HDAC7 around the non-phosphorylation site of S155.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Isotype</b>	IgG
<b>Alias</b>	HDAC7; HDAC7A; Histone deacetylase 7; HD7; Histone deacetylase 7A; HD7a
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	HDAC7

## Image



Western Blot analysis of PC-3 cells using HDAC7 Polyclonal Antibody



Western Blot analysis of PC3 KB MCF7 HELA cells using HDAC7 Polyclonal Antibody

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.