




HRH1 Antibody

Product Code	CSB-PA002911
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35367
Immunogen	Synthesized peptide derived from the Internal region of Human Histamine H1 Receptor.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	HRH1; Histamine H1 receptor; H1R; HH1R
Immunogen Species	Homo sapiens (Human)
Target Names	HRH1

Image



Western Blot analysis of COLO205 cells using Histamine H1 Receptor Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.