




HDAC7A Antibody

Product Code	CSB-PA002916
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WUI4
Immunogen	Synthesized peptide derived from the C-terminal region of Human Histone deacetylase 7a.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	HDAC7; HDAC7A; Histone deacetylase 7; HD7; Histone deacetylase 7A; HD7a
Immunogen Species	Homo sapiens (Human)
Target Names	HDAC7

Image



Western Blot analysis of HuvEc cells using Histone deacetylase 7a Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.